


Derma Free[®]

POWDER-FREE VINYL EXAMINATION GLOVES

• STANDARD EXAMINATION GLOVE

Backed by the industry's leading quality system, Derma Free[®] powder-free vinyl is 100% synthetic and features a smooth surface for enhanced tactile sensitivity.

Not made from natural rubber latex.



Comfort, performance and value in a powder-free medical examination glove.

Backed by the industry's leading quality system, **Derma Free[®]** powder-free vinyl gloves provide solid hand protection for a variety of applications. Thicker than typical vinyl gloves, the **Derma Free[®]** glove is 100% synthetic and features a smooth surface for more natural tactile sensitivity.

MICROFLEX[®]
THE MOST TRUSTED NAME IN GLOVES[®]

REGISTERED
ISO13485

Derma Free[®]

POWDER-FREE VINYL EXAMINATION GLOVES

STANDARD EXAMINATION GLOVE


GLOVE PROFILE

	ASTM*	Derma Free [®]
Length (mm/in):	230	230/9.05
Cuff Thickness (mils):	N/A	3.5
Palm Thickness (mils):	min. 3.1	4.3
Finger Thickness (mils):	min. 2.0	3.9
Tensile Strength		
Before Aging (MPa):	min. 11	min. 11
After Aging (MPa):	min. 11	min. 11
Elasticity/Elongation (%):		
Before Aging (MPa):	min. 300	min. 350
After Aging (MPa):	min. 300	min. 300

GLOVE ATTRIBUTES

Material:	Vinyl; not made from natural rubber latex
Type:	Non-sterile
Shape:	Ambidextrous
Color:	Clear
Sizes:	X-Small, Small, Medium, Large, X-Large
Cuff:	Beaded
Glove Interior:	Powder-free
Glove Exterior:	Smooth
Usage:	For single use only
Packaging:	100 gloves per dispenser box, by weight 10 Boxes per case

*ASTM D 5250



Powder-free Vinyl Exam Gloves

See for yourself why **Microflex[®]** gloves are the first choice of examination glove users nationwide. For samples or more information on any of the fine products from **Microflex[®]**, contact your **Microflex[®]** distributor. For the name of a distributor in your area, call:

(800) 876-6866

Or visit our website at:

www.microflex.com



PRODUCT CODE	RANGE OF SIZE	DESCRIPTION	CASE QUANTITY
DF-850-XS	EXTRA-SMALL	Vinyl, Powder-free	1000
DF-850-S	SMALL	Vinyl, Powder-free	1000
DF-850-M	MEDIUM	Vinyl, Powder-free	1000
DF-850-L	LARGE	Vinyl, Powder-free	1000
DF-850-XL	EXTRA-LARGE	Vinyl, Powder-free	1000

Caution Statement: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.