

KIMBERLY-CLARK\* **SAFESKIN**\*  
 PURPLE NITRILE\*,  
 PURPLE NITRILE-XTRA\*,  
 PURPLE NITRILE\*  
 Sterile Singles & Pairs Exam Gloves



As a healthcare professional, you know the importance of protecting yourself from daily exposure to bacteria, viruses and chemicals. Just as the proper glove creates a barrier between you and these dangerous elements, it also enhances your ability to offer patients maximum security against cross-contamination.

From proprietary latex and nitrile formulations for moderate to high risk procedures involving infectious fluids or bloodborne pathogens, to vinyl for low risk situations, KIMBERLY-CLARK\* **SAFESKIN**\* has a glove solution for you. Our complete line includes gloves for primary care services, dental procedures, diagnostic tests, surgical procedures, and general use applications.

**Functional Benefits:**

- Unique purple color provides quick, visual differentiation from natural rubber latex and other synthetic gloves.
- Textured fingertips enhance wet and dry grip.
- Custom design enhances comfort and fit.
- Available in 9½" and 12" lengths.
- Available in sterile singles and pairs.

**Clinical Specifications:**

- Natural Rubber Latex-free – primary material is acrylonitrile-butadiene, reducing the potential for glove-associated Type I allergic reactions to natural rubber latex proteins.
- Powder-free, minimizing the potential for powder-related complications, such as irritant contact dermatitis.
- PURPLE NITRILE-XTRA® material has been tested for use with ten chemotherapy drugs and gluteraldehyde to determine permeation and breakthrough times. Test results on file.

## Kimberly-Clark\* SAFESKIN® Exam Gloves

### Product Specifications:

PURPLE NITRILE® Measurements (Target)	gauge thickness (measuring single thickness)		average length
	MM	MIL	
Middle Finger (15mm from tip)	0.15	5.9	242 mm or 9.5" 305 mm or 12" <i>(measured from the tip of the middle finger to the edge of the cuff)</i>
Palm	0.12	4.7	
Cuff	0.09	3.5	

### Physical Properties (Target)

before aging		after aging	
Tensile Strength	Ultimate Elongation	Tensile Strength	Ultimate Elongation
21 MPa	550%	21 MPa	500%

Note: all specifications are subject to change without notice.

KIMBERLY-CLARK\* SAFESKIN®  
PURPLE NITRILE\*,  
PURPLE NITRILE-Xtra\*,  
PURPLE NITRILE\*  
Sterile Singles & Pairs  
Exam Gloves



### Quality Standards:

- Exceeds current ASTM D 6319 standard for critical defects (AQL-2.5). Safeskin AQL for critical defects is 1.5.
- Manufactured in accordance with Quality System ISO 9002.
- Biocompatible as measured by Primary Skin Irritation, Repeat Challenge Sensitization and 200 Person Modified Draize Test.
- Sterilized in accordance with AAMI guidelines at a Sterility Assurance Level (SAL) of 10<sup>-6</sup>.

From protective apparel to medical devices – everything we make, everything we do, is focused on protecting healthcare professionals and patients. This means we work every day to develop and deliver quality products and to provide responsive customer service that you can depend on to meet your protection needs.

### KIMBERLY-CLARK\* SAFESKIN® PURPLE NITRILE® – 9.5"

Stock #	Size	Gloves/ Box	Boxes/ Case
55080	X-Small	100	10
55081	Small	100	10
55082	Medium	100	10
55083	Large	100	10
55084	X-Large	90	10

### KIMBERLY-CLARK\* SAFESKIN® PURPLE NITRILE-XTRA® – 12"

Stock #	Size	Gloves/ Box	Boxes/ Case
55090	X-Small	50	10
50601	Small	50	10
50602	Medium	50	10
50603	Large	50	10
50604	X-Large	50	10

### PURPLE NITRILE® – Sterile Singles

Singles Stock #	Size	Gloves/ Box	Boxes/ Case
52101	Small	100	4
52102	Medium	100	4
52103	Large	100	4

### PURPLE NITRILE® – Sterile Pairs

Stock #	Size	Pairs/ Box	Boxes/ Case
55091	Small	50	4
55092	Medium	50	4
55093	Large	50	4



### Commitment to Excellence

If, for any reason, our products do not meet your expectations, please let us know your comments or suggestions for improvement. Your input will result in a concerted effort on our part to meet your requirements. Our goal is to provide quality products that completely meet your needs time after time.

For more information, please call 1-800-KCHELPS in the United States, or visit our web site at [www.kchealthcare.com](http://www.kchealthcare.com).

\* Registered Trademark or Trademark of Kimberly-Clark Corp. Roswell, GA 30076 or its affiliates.  
© The COLOR PURPLE is a Registered Trademark of Kimberly-Clark Corporation.  
Manufactured for Kimberly-Clark Corp. by SAFESKIN Group (Thailand), a wholly owned subsidiary of Kimberly-Clark Corporation, Roswell, GA USA. © 2004 KCWW. All rights reserved.

Printed on CLASSIC CREST® Text, Avon Brilliant White, 100 lb. produced by the Neenah Paper Division of Kimberly-Clark Corporation.

KLM-92-1 H7174

Contact your  
Kimberly-Clark representative  
to schedule a glove evaluation  
or to place an order.

Protection. For life.

