


Kimberly-Clark® **SAFESKIN®** Medical Gloves

KIMBERLY-CLARK® SAFESKIN®

Synthetic Plus Vinyl Exam Gloves



As a healthcare professional, you know the importance of protecting yourself from daily exposure to bacteria, viruses and chemicals. Just as the proper glove creates a barrier between you and these dangerous elements, it also enhances your ability to offer patients maximum security against cross-contamination.

From proprietary latex and nitrile formulations for moderate to high risk procedures involving infectious fluids or bloodborne pathogens, to vinyl for low risk situations, KIMBERLY-CLARK® **SAFESKIN®** has a glove solution for you. Our complete line includes gloves for primary care services, dental procedures, diagnostic tests, surgical procedures, and general use applications.

Functional Benefits:

- Preferred fit and feel compared to SAFESKIN® Synthetic clear vinyl gloves.
- White color provides a visual cue that the glove is made from synthetic materials.
- Polymer coating improves donning.
- Ambidextrous for use across a variety of glove applications.

Clinical Specifications:

- Natural Rubber Latex-free – primary material is polymer coated polyvinyl chloride, reducing the potential for glove-associated Type I allergic reactions to natural rubber latex proteins.
- Powder-free, minimizing the potential for powder-related complications, such as irritant contact dermatitis.

Protection. For life.

 **Kimberly-Clark**

Kimberly-Clark® **SAFESKIN® Medical Gloves**

KIMBERLY-CLARK® SAFESKIN®

Synthetic Plus Vinyl Exam Gloves



Product Specifications:

<i>Synthetic Plus Measurements (Target)</i>	gauge thickness <i>(measuring single thickness)</i>		average length <i>(measured from the tip of the middle finger to the edge of the cuff)</i>
	MM	MIL	
Middle Finger <i>(15mm from tip)</i>	0.20	7.9	245 mm or 9.5"
Palm	0.13	5.1	
Cuff	0.11	4.3	

Physical Properties (Target)

before aging	
<i>Tensile Strength</i>	<i>Ultimate Elongation</i>
12 MPa	350%

Note: all specifications are subject to change without notice

Quality Standards:

- Exceeds current ASTM D 5250 standard for critical defects (AQL-2.5). Safeskin AQL for critical defects is 1.5.

From protective apparel to medical devices – everything we make, everything we do, is focused on protecting healthcare professionals and patients. This means we work every day to develop and deliver quality products and to provide responsive customer service that you can depend on to meet your protection needs.

KIMBERLY-CLARK® SAFESKIN® Synthetic Plus Vinyl Exam Gloves

Stock #	Size	Gloves/Box	Boxes/Case
50031	Small	100	10
50032	Medium	100	10
50033	Large	100	10
50034	X-Large	90	10



Commitment to Excellence

If, for any reason, our products do not meet your expectations, please let us know your comments or suggestions for improvement. Your input will result in a concerted effort on our part to meet your requirements. Our goal is to provide quality products that completely meet your needs time after time.

For more information, please call 1-800-KCHELPS in the United States, or visit our web site at www.kchealthcare.com.

* Registered Trademark or Trademarks of Kimberly-Clark Corp. Roswell, GA 30076 or its affiliates. Manufactured for Kimberly-Clark Corp. by SAFESKIN Group (Thailand), a wholly owned subsidiary of Kimberly-Clark Corporation, Roswell, GA USA. © 2001 KCC. All rights reserved.

Printed on CLASSIC CREST® Text, Avon Brilliant White, 100 lb. produced by Kimberly-Clark's Neenah Paper.

KLM 22 H6100