

Aloetouch® Ultra IC



Someday all exam gloves
will be made this way

Thanks to Medline's new proprietary inner surface coating, our Aloetouch® Ultra IC latex free exam glove is faster and easier to don,* even when your hands are wet.

The new inner surface coating combines the moisturizing benefits of pure aloe vera gel with an advanced polymer for trouble-free donning.

Ask your Medline rep to show you this extraordinary glove.



*Independent Medline study.
Data available upon request.

NEW!

Aloetouch® Ultra IC


Powder-Free Stretch Synthetic Coated Exam Gloves

To make the already popular Aloetouch® Ultra exam glove even better, industry leader Medline went to the people who know gloves best—nurses!

The result—the Aloetouch® Ultra IC retains our patented freeze-dried aloe lining and couples it with a proprietary enhanced inner coating. Pairing these features means the Aloetouch® Ultra IC doesn't just moisturize the hands with soothing aloe vera—it's faster and easier to don.*



Enhanced Inner Coating



Soothing Aloe Vera



latexfree
PRODUCT

Aloetouch® Ultra IC High-Tech Features

- Proprietary inner coating makes Aloetouch® Ultra IC gloves easier and faster to don
- The pure, freeze-dried aloe vera in the coating softens hands with each use
- 100 percent latex-free protection against bloodborne pathogens
- DEHP-and accelerator-free
- Easy-on tapered cuffs for a better fit

To guarantee Aloetouch® Ultra IC exam gloves stand up to nurses' demands, Medline put them on their hands! Field trials were conducted at hospitals currently using the Aloetouch® Ultra exam gloves.*

Nurses said the new Aloetouch® Ultra IC...

"Is easier to don when my hands are wet."

"Is stronger and more resistant to tearing."

"Offers enhanced tactile sensitivity."

"Is noticeably better than the gloves I normally wear."

"Fits my hand more comfortably."

The Aloetouch® Ultra IC—someday, all exam gloves will be made this way.

Item Number	Size	Packaging
MDS195073	X-Small	100 ea/bx, 10bx/cs
MDS195074	Small	100 ea/bx, 10bx/cs
MDS195075	Medium	100 ea/bx, 10bx/cs
MDS195076	Large	100 ea/bx, 10bx/cs
MDS195077	X-large	100 ea/bx, 10bx/cs

Specification	ASTM D5250 (Min.)	Aloetouch Ultra IC (Min.)
Color/Texture	N/A	Green/Smooth
Glove Length	230mm	240mm
Palm Thickness	0.08mm	0.12mm
Finger Thickness	0.05mm	0.14mm
Tensile Strength Before Aging	11 MPa	13 MPa
Elongation Before Aging	300%	500%
Viral Penetration Test	N/A	PASSED
AQL	FDA 4.0	2.5

*Independent Medline study. Data available upon request.

©2007 Medline Industries, Inc. Medline and Aloetouch are registered trademarks of Medline Industries, Inc. We reserve the right to correct any errors that may occur within this sell sheet.

MKT207003/LIT542/5M/WMP8



Medline Industries, Inc.
One Medline Place
Mundelein, IL 60060

1-800-MEDLINE
www.medline.com
www.aloetouch.com

Medline Canada
1-800-396-6996
www.medline.ca
canada@medline.com

Not all products are currently available for sale in Canada. Please contact your account representative for additional information.