



Sensicare[®] Nitrile

- **Strong**
- **Comfortable**
- **Affordable**
- **Chemo Tested**

Revolutionizing the nitrile marketplace with greater elasticity and exceptional barrier protection



Sensicare Nitrile Powder-Free Exam Gloves

Medline has leveraged advanced glove technology to produce a new Powder-Free Nitrile glove that will revolutionize the way that nitrile gloves are viewed in the healthcare marketplace. In the past, nitrile gloves were thick and rigid and typically only used in high-risk areas. The new Sensicare breaks this mold as a soft yet strong nitrile glove that provides exceptional barrier protection – it has been tested for use with chemotherapy drugs!

The Sensicare is an excellent option for Healthcare facilities that are looking to move towards a latex-free environment. Healthcare professionals now have a nitrile option that can be easily used facility-wide. Clinicians marvel at the outstanding tactile sensitivity and excellent fit that Sensicare provides.



Features

100% Latex-Free

Strong and Durable

Soft and supple for excellent dexterity

High elasticity

Tested for use with chemotherapy drugs *

Fingertip Textured

Great feel for I.V. starts

Fresh, blue color

Powder-Free

Item Number	Size	Packaging
484800	X-Small	100 ea/bx, 10bx/cs
484801	Small	100 ea/bx, 10bx/cs
484802	Medium	100 ea/bx, 10bx/cs
484803	Large	100 ea/bx, 10bx/cs
484804	X-Large	100 ea/bx, 10bx/cs

Specification	ASTM*	Medline
Viral Penetration Test	N/A	Passed
Palm Thickness (mm)	.05 min	.09mm
Finger Thickness (mm)	.05 mm	.11mm
Tensile Strength Before Aging	12.5 Mpa	26.6 Mpa
Elongation Before Aging	500%	531%
Color/Texture	N/A	Blue/Fingertip Textured



* Please contact Medline for questions regarding chemotherapy drug testing

ASTM Standard D6319

1-800-MEDLINE (633-5463) / www.medline.com

©2005 Medline Industries, Inc. One Medline Place, Mundelein, IL 60060 Medline and Sensicare are registered trademarks of Medline Industries, Inc
MKT205252/LIT420/10M/4WMP