

A close-up photograph of a person's hands wearing white, latex-free exam gloves. The person is wearing a yellow surgical gown. The background is a solid blue color.

“Finally!

We've spent months looking for a glove that allows us to go latex-free. Medline's Universal 3G is the first glove that gives the comfort, barrier protection and price we need to make the switch. We're converting our entire health system to the Universal 3G. //

Desiree Roderick
Sr. Purchasing Manager
Valley Health System
Hemet, CA

Medline's
UNIVERSAL 3G™
Exam Glove

Introducing...

Universal 3G™ Exam Gloves

Only from Medline

Medline has advanced the technology of stretch synthetic exam gloves. The result is a patented glove formulation that delivers a softness no one thought possible, with outstanding strength and durability. The Universal 3G is the first synthetic glove to offer truly latex-like performance.

The resin in the Universal 3G glove has much longer polymer chains than other synthetic gloves on the market today. The longer polymer chains make the resin MORE workable, MORE flexible and MORE resilient. Now you can standardize on a single latex-free glove.



Made from a proprietary formula, the Universal 3G provides the tactile sensitivity, fit and barrier protection that gives facilities across the country the freedom to go latex-free!

Features

100% Latex-Free

Soft and supple for excellent dexterity

Great "feel" for I.V. starts

Odor free and accelerator free

Tapered cuff for a better fit

Clean, white color

Powder-Free

Non-sterile

Exceeds ASTM D5250

Item Number	Size	Packaging
MDS196073	X-Small	100 ea/bx, 10bx/cs
MDS196074	Small	100 ea/bx, 10bx/cs
MDS196075	Medium	100 ea/bx, 10bx/cs
MDS196076	Large	100 ea/bx, 10bx/cs
MDS196077	X-Large	100 ea/bx, 10bx/cs

Specification	ASTM*	Universal 3G**
Color/Texture	N/A	White
Glove Length	230mm	240mm
Palm Thickness (mm)	.08mm	.13mm
Finger Thickness (mm)	.08mm	.15mm
Tensile Strength Before Aging	9 Mpa	14.8 Mpa
Elongation Before Aging	300%	520%
Viral Penetration Test	N/A	Passed

*ASTM D5250-00

** Mean values. Test results on file

Patent Pending



1-800-MEDLINE (633-5463) / www.medline.com

©2004 Medline Industries, Inc. Medline is a registered trademark and Universal 3G is a trademark of Medline Industries, Inc. One Medline Place, Mundelein, IL 60060

MKT204011/LIT136/20M