

# MICRO ONE®

LIGHTLY POWDERED LATEX EXAMINATION GLOVES

• STANDARD EXAMINATION GLOVE



## Features & Benefits:

- Longer than typical latex examination gloves for added protection where you need it most
- Lightly powdered with absorbable cornstarch (USP) for easy donning
- Smooth surface for excellent tactile sensitivity
- Meets federal guidelines for maximum powder and maximum water extractable protein content



## Lightly powdered. A Proven Performer.

The **Micro One**® glove is a highly versatile examination glove that performs well in demanding conditions. Made from soft latex, **Micro One**® gloves have a natural fit for reduced hand fatigue. Lightly powdered **Micro One**® gloves don easily and have a smooth surface for uncompromised sensitivity. **Micro One**® gloves are thicker and longer than typical latex exam gloves for added protection where you need it most. The **Micro One**® glove meets or exceeds current ASTM standards for exam glove quality and performance.

**MICRO FLEX**®  
THE MOST TRUSTED NAME IN GLOVES®

REGISTERED  
5013485

# MICRO ONE®

LIGHTLY POWDERED LATEX EXAMINATION GLOVES

● STANDARD EXAMINATION GLOVE



## GLOVE PROFILE

	ASTM*	Micro One®
<b>Length (mm/in):</b>	220 / 8.7	245 / 9.4
<b>Cuff Thickness (mils):</b>	n/a	3.9
<b>Palm Thickness (mils):</b>	min. 3.1	5.5
<b>Finger Thickness (mils):</b>	min. 3.1	6.3
<b>Tensile Strength</b>		
Before Aging (MPa):	min. 18	min. 21
After Aging (MPa):	min. 14	min. 18
<b>Elasticity/Elongation:</b>		
Before Aging (%):	min. 650	750
After Aging (%):	min. 500	600
<b>Reduced Protein Rating:</b>	max. 200 µg/dm <sup>2</sup>	max. 150 µg/g

## GLOVE ATTRIBUTES

<b>Material:</b>	Premium quality natural rubber latex
<b>Shape:</b>	Ambidextrous
<b>Color:</b>	Natural
<b>Sizes:</b>	X-Small, Small, Medium, Large, X-Large
<b>Cuff:</b>	Beaded
<b>Glove Interior:</b>	Lightly powdered with absorbable cornstarch, USP
<b>Glove Exterior:</b>	Smooth, non-textured
<b>Usage:</b>	For single use only
<b>Packaging:</b>	100 gloves per dispenser box, by weight 10 Boxes per case

\* ASTM D 3578

## Latex Examination Gloves

See for yourself why **Microflex®** gloves are the first choice of exam glove users nationwide. For samples or more information on any of the fine latex and synthetic products from **Microflex®**, contact your **Microflex®** distributor or for the name of a distributor in your area, call:

**(800) 876-6866**

Or visit our website at:

[www.microflex.com](http://www.microflex.com)

Fax: (800) 876-6632.



PRODUCT CODE	RANGE OF SIZE	DESCRIPTION	CASE QUANTITY
MO-150-XS	EXTRA-SMALL	Latex, Lightly Powdered, Examination	1000
MO-150-S	SMALL	Latex, Lightly Powdered, Examination	1000
MO-150-M	MEDIUM	Latex, Lightly Powdered, Examination	1000
MO-150-L	LARGE	Latex, Lightly Powdered, Examination	1000
MO-150-XL	EXTRA-LARGE	Latex, Lightly Powdered, Examination	1000

**CAUTION STATEMENT:** This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

**MICRO FLEX®**  
THE MOST TRUSTED NAME IN GLOVES®  
2301 ROBB DRIVE • RENO, NV 89523