



**BLACK POWDER-FREE NITRILE
STANDARD EXAMINATION GLOVES**

Recommended for use in the following industries:

Laboratory / Law Enforcement / Industrial Safety / Tattoo /
Automotive / Corectional Facilities / Funeral

When enhanced tactile sensitivity and consistent grip are a must, make MidKnight™ your glove of choice.

MidKnight™ gloves offer the superior combination of durability and protection against day-to-day hazards, in a comfortable nitrile glove. The glove of choice for individuals who demand uncompromising quality and protection, MidKnight™ provides an exceptional balance of performance, comfort and durability. Designed specifically with a distinct color that resists stains.

Recommended for use when:

- Day-to-day hazards require the best combination of durability and consistency
- Working in environments that demand firm grip and reliability
- Precise contact requires enhanced tactile sensitivity

NOT MADE FROM NATURAL RUBBER LATEX

MidKnight™ Attributes



For current specification information, please visit www.microflex.com.

| | ASTM D 6319 | MidKnight™ |
|---------------------------------|-------------|------------|
| Length (mm/in): | 220 / 8.7 | 245 / 9.6 |
| Cuff Thickness (mils): | n/a | 3.1 |
| Palm Thickness (mils): | min. 1.9 | 4.7 |
| Finger Thickness (mils): | min. 1.9 | 5.5 |
| Tensile Strength | | |
| Before Aging (MPa): | min. 14 | 16 |
| After Aging (MPa): | min. 14 | 14 |
| Elasticity/Elongation: | | |
| Before Aging (%): | min. 500 | 550 |
| After Aging (%): | min. 400 | 400 |

Order #

XS... MK-296-XS L... MK-296-L
S... MK-296-S XL... MK-296-XL
M... MK-296-M XXL... MK-296-XXL

Material: Nitrile
Type: Non-sterile
Color: Black
Cuff: Beaded
Glove Exterior: Fully Textured (XXL glove is fingertip textured only; not fully textured)
Packaging: 1000 gloves per case, by weight, 10 Boxes per case
100 gloves per dispenser box, by weight

Usage: For single use only
Shape: Ambidextrous
Size: XS - XXL
Glove Interior: Powder-free

Caution: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.
U.S. Patent No. Re. 35,616, European Patent No. 0456333

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. Contact a Microflex® distributor today to purchase quality Microflex® products. For the name of a distributor in your area, call: **(800) 876-6866**. Fax: (800) 876-6632.

