



Showa Best Glove will donate a portion of the proceeds from N-DEX BWC to organizations promoting breast cancer research and awareness.

NITRILE MEDICAL EXAM GLOVES



BECAUSE WE CARE



6205PF glove

Disposable Exam Glove

Characteristics

- 100% Nitrile Glove
- Powder Free
- Ambidextrous
- 4 mils, 9.5" Long
- Textured Fingertips
- 10 Dispensers/100 gloves ea.

Sizes:

6205PFXS	Extra Small
6205PFS	Small
6205PFM	Medium
6205PFL	Large
6205PFXL	Extra Large

Best N-DEX products are the only disposable nitrile gloves manufactured in the USA.

Description

Showa Best Glove's N-DEX® BWC™ (Because We Care) PINK® nitrile glove enables disposable glove wearers to show their support of cancer survivors while enjoying all the characteristics associated with the Best N-DEX hand protection. These 100% nitrile, non-latex medical grade gloves are part of the world renowned N-DEX family of disposable hand protection proudly manufactured in America. N-DEX BWC gloves are ambidextrous, 9.5" long, 4-mils thick with textured fingertips and are packed 10 dispensers of 100 gloves per case.

Applications/Certifications

- Blood-Borne Pathogen Protection
- Patient Exam
- Chemotherapy Administration
- Laboratories
- Class One Medical Glove
- NFPA 1999 Emergency Medical Operations
- In Compliance 21 CFR, Parts 170-199

This list of suggested applications is partial and based on laboratory testing of glove performance.

For a complete list of suggested applications, visit our website:

www.showabestglove.com

PINK® is a registered trademark of Colour World, LLC.