

SHOWA 

Best 

N-DEX[®] Free

**N-DEX[®] Free Ultimate
9500PF**

**100% NITRILE GLOVE
ACCELERATOR-FREE
DISPOSABLE 6 MIL, 11" GLOVE
NFPA 1999 COMPLIANT EMS GLOVE
CLASS ONE MEDICAL GLOVE
US PATENT 35.616**



Disposable Glove:

Best Part Number	Size	Packaging
9500PFXS	Extra-Small	50 Gloves per Dispenser 20 Dispensers per Case
9500PFS	Small	
9500PFM	Medium	
9500PFL	Large	
9500PFXL	Extra-Large	

Description

SHOWA Best Inc's N-DEX[®] Free Ultimate 9500PF is longer and thicker. The Accelerator-Free Disposable 100% Nitrile gloves offers the Ultimate in dexterity, comfort and NFPA 1999 Certified protection. For every-day tasks to life-saving Emergency Medical Operations, Showa Best's Ultimate N-DEX Free disposable provides the protection needed on the hands of EMT, painters, automotive and household workers with barrier protection, chemical and puncture resistance and tactile discrimination. The 5 mil, 11 inch long disposable nitrile glove provides wearers barrier protection and comfort unequalled by supported nitrile gloves nearly twice its thickness. Ultimate N-DEX Free glove's patented construction includes Best's low-modulus formulation, maximizing comfort and minimizing hand fatigue.

Cat. III

EN 388

0000

Performance Levels:

Abrasion: Level 0-4
Cut Resistance: Level 0-5
Tear: Level 0-4
Puncture: Level 0-4

EN-374-3 Chemical Resistance



Diesel Fuel 0 1 2 3 4 5 **6**
Sodium Hydroxide 40% 0 1 2 3 4 5 **6**
Hydrochloric Acid 37% 0 1 2 3 4 5 **6**

ANSI/ISEA 105-2005 Performance

Abrasion 0 1 **2** 3 4 5 6
Cut Resistance 0 **1** 2 3 4 5
Puncture 0 **1** 2 3 4 5

Suggested Applications

- Emergency Medical Operations
- Manufacturing
- Food Processing
- Laboratories
- Food Handling
- Janitorial
- Paint Applications
- Fuel Handling
- Chemotherapy Administration

This list of suggested applications is partial and based on laboratory testing of glove performance. For a complete list of suggested applications, visit our website:

www.showabestglove.com

USA:
Best Glove
579 Edison Street
Menlo, GA 30731
Phone: 800-241-0323
Fax: 888-393-2666
Email: USA@bestglove.com

Canada:
Best Glove Manufacturing, Ltd.
253 Michaud Street
Coaticook, Quebec J1A 1A9 Canada
Phone: 800-565-BEST
Fax: 819-849-6120
Email: Canada@bestglove.com

Europe:
Best Manufacturing Europe N.V.
Draaiboornstraat 6 bus 12
2160 Wommelgem Belgium
Phone: +32 (0)3 458 33 33
Fax: +32 (0)3 458 34 57
Email: Europe@bestglove.com

Australasia:
Best Manufacturing Australia Pty. Ltd.
PO Box 187, Bexley
NSW 2207 Australia
Phone: +61 (0)2 8509 5401
Fax: +61 (0)2 8509 5403
Email: Australiansales@bestglove.com