



N-DEX[®] Free

N-DEX[®] Free 7705PFT



- 100% NITRILE GLOVE
- ACCELERATOR-FREE GLOVE
- DISPOSABLE 4 MIL, 9.5" GLOVE
- NFPA 1999 COMPLIANT EMS GLOVE
- CLASS ONE MEDICAL GLOVE
- US PATENT 35.616
- US PATENT 6706816

Disposable Glove:

Part Number	Size	Packaging
7705PFTXS	Extra-Small	100 Gloves per Dispenser 20 Dispensers per Case
7705PFTS	Small	
7705PFTM	Medium	
7705PFTL	Large	
7705PFTXL	Extra-Large	

Description

Best's revolutionary N-DEX[®] Free 7705PFT is the world's first accelerator free disposable nitrile glove. Free from the accelerators common in all disposable nitrile gloves, which have been linked to workplace conditions like Type IV dermatitis. Free from Natural Rubber Latex and Type I allergic reactions linked to Latex proteins. Developed by Best to the same exacting standards of its N-DEX family of hand protection, N-DEX Free's patented accelerator free formulation creates a luminous, fun, yet PPE-conscious, high visibility green glove. The hi-vis color is great for jobs where increased visibility and color-coded operations are critical. Available with textured fingertips, N-DEX Free provides laboratory, pharmaceutical and industrial workers an accelerator-free, Latex-free option with greater grip, comfort and dexterity.

Suggested Applications

- Emergency Medical Operations
- Manufacturing
- Food Processing
- Laboratories
- Janitorial
- Paint Applications
- Fuel Handling
- Patient Examination

This list of suggested applications is partial and based on laboratory testing of glove performance. For a complete list of suggested applications, visit our website:

www.showabestglove.com

EN 455 EN 374-3



AQL of 1.5

EN-374-3 Chemical Resistance



Diesel Fuel 0 1 2 **3** 4 5 6
Sodium Hydroxide 40% 0 1 2 3 4 **5** 6
Hydrochloric Acid 37% 0 1 2 3 **4** 5 6

USA:
Showa Best Glove
579 Edison Street
Menlo, GA 30731
Phone: 800-241-0323
Fax: 888-393-2666
Email: USA@showabestglove.com

Canada:
Showa Best Glove Manufacturing, Ltd.
253 Michaud Street
Coaticook, Quebec J1A 1A9 Canada
Phone: 800-565-BEST
Fax: 819-849-6120
Email: Canada@showabestglove.com

Europe:
Showa Best Glove
Tour Franklin
La Defense 8
92042 Paris La Defense Cedex France
Tel Direct +33 (0)1 55 62 14 20
Fax + 33 (0)1 55 62 14 29
Email: Info@showabestglove.eu

Australasia:
Best Manufacturing Australia Pty. Ltd.
PO Box 187, Bexley
NSW 2207 Australia
Phone: +61 (0)2 8509 5401
Fax: +61 (0)2 8509 5403
Email: Australiansales@bestglove.com