



# N-DEX® # 6005PF

100% NITRILE GLOVE  
DISPOSABLE 4 MIL, 9.5" GLOVE  
NFPA 1999 COMPLIANT EMS GLOVE  
CLASS ONE MEDICAL GLOVE

## Disposable Glove:

Part Number	Size	Packaging
6005PFXS	Extra-Small	100 Gloves per Dispenser 20 Dispensers per Case
6005PFS	Small	
6005PFM	Medium	
6005PFL	Large	
6005PFXL	Extra-Large	

## Description

**SHOWA Best Inc's N-DEX® 6005PF glove** is a 100% Nitrile, non-Latex, powder-free medical grade glove which features Best's proprietary beaded and rolled cuff design. The N-DEX® 6005PF are ambidextrous, 9 1/2" long, 4-mil thick, packed 2,000 gloves per case. The Class One Medical Examination glove meets OSHA guidelines for handling blood-borne pathogens. The glove is acceptable for food handling applications and is in compliance with 21CFR, parts 170-199. N-DEX 6005PF gloves are manufactured in an FDA-inspected and ISO 9002 registered facility. N-DEX 6005PF gloves are certified as compliant with NFPA 1999 Standard for Emergency Medical Operations and they are made in the USA.

## Suggested Applications

- Emergency Medical Operations
- Manufacturing
- Food Processing
- Laboratories
- Janitorial
- Paint Applications
- Fuel Handling
- Patient Examination

This list of suggested applications is partial and based on laboratory testing of glove performance. For a complete list of suggested applications, visit our website:

[www.showabestglove.com](http://www.showabestglove.com)

EN 455 EN 374-3



AQL of 1.5



### EN-374-3 Chemical Resistance

Diesel Fuel	Sodium Hydroxide 40%
0 1 2 3 4 5 6	0 1 2 3 4 5 6
Hydrochloric Acid 37%	
0 1 2 3 4 5 6	

USA:  
Best Glove  
579 Edison Street  
Menlo, GA 30731  
Phone: 800-241-0323  
Fax: 888-393-2666  
Email: [USA@bestglove.com](mailto:USA@bestglove.com)

Canada:  
Best Glove Manufacturing, Ltd.  
253 Michaud Street  
Coaticook, Quebec J1A 1A9 Canada  
Phone: 800-565-BEST  
Fax: 819-849-6120  
Email: [Canada@bestglove.com](mailto:Canada@bestglove.com)

Europe:  
Best Manufacturing Europe N.V.  
Draaiboornstraat 6 bus 12  
2160 Wommelgem Belgium  
Phone: +32 (0)3 458 33 33  
Fax: +32 (0)3 458 34 57  
Email: [Europe@bestglove.com](mailto:Europe@bestglove.com)

Australasia:  
Best Manufacturing Australia Pty. Ltd.  
PO Box 187, Bexley  
NSW 2207 Australia  
Phone: +61 (0)2 8509 5401  
Fax: +61 (0)2 8509 5403  
Email: [Australiansales@bestglove.com](mailto:Australiansales@bestglove.com)