



# Ultimate N-DEX<sup>®</sup> B9905PF

## Clean-Processed Disposable Nitrile Gloves

**CLEAN-PROCESSED  
100% NITRILE GLOVE  
DISPOSABLE, 6 MIL, 11" GLOVE  
NFPA 1999 COMPLIANT EMS GLOVE  
CLASS ONE MEDICAL GLOVE  
US PATENT 35,616**



### Disposable Glove:

Best Part Number	Size	Packaging
B9905PFXS	X-Small	50 Gloves per Polybag  20 Polybags per Case
B9905PFS	Small	
B9905PFM	Medium	
B9905PFL	Large	
B9905PFXL	Extra-Large	

### Description

Best's Clean-Processed N-DEX<sup>®</sup> Ultimate longer, thicker disposable 100% Nitrile gloves offer the Ultimate in dexterity, comfort and NFPA 1999 Certified protection. Packed in Cleanroom double bags, these products may be suitable for Class 1,000 Cleanrooms. Clean-Processed Ultimate N-DEX disposable provides the protection needed on the hands of EMT, painters, automotive and Cleanroom personnel with barrier protection, chemical and puncture resistance and tactile discrimination. The 6 mil, 11" long disposable nitrile glove provides wearers barrier protection and comfort unequalled by supported nitrile gloves nearly twice its thickness. Ultimate N-DEX glove's patented construction includes Best's low-modulus formulation, maximizing comfort and minimizing hand fatigue.

### Suggested Applications

- Cleanroom Operations
- Manufacturing
- Food Processing
- Laboratories
- Electronics
- Paint Applications
- Fuel Handling
- Chemotherapy Administration

This list of suggested applications is partial and based on laboratory testing of glove performance. For a complete list of suggested applications, visit our website:

[www.showabestglove.com](http://www.showabestglove.com)

EN 455



### ANSI/ISEA 105-2005 Performance

**Abrasion** 0 1 **2** 3 4 5 6      **Cut Resistance** 0 **1** 2 3 4 5      **Puncture** 0 **1** 2 3 4 5

USA:  
Showa Best Glove  
579 Edison Street  
Menlo, GA 30731  
Phone: 800-241-0323  
Fax: 888-393-2666  
Email: [USA@showabestglove.com](mailto:USA@showabestglove.com)

Canada:  
Showa Best Glove Manufacturing, Ltd.  
253 Michaud Street  
Coaticook, Quebec J1A 1A9 Canada  
Phone: 800-565-BEST  
Fax: 819-849-6120  
Email: [Canada@showabestglove.com](mailto:Canada@showabestglove.com)

Europe:  
Showa Best Glove  
Tour Franklin  
La Defense 8  
92042 Paris La Defense Cedex France  
Tel Direct +33 (0)1 55 62 14 20  
Fax + 33 (0)1 55 62 14 29  
Email: [Info@showabestglove.eu](mailto:Info@showabestglove.eu)

Australasia:  
Best Manufacturing Australia Pty. Ltd.  
PO Box 187, Bexley  
NSW 2207 Australia  
Phone: +61 (0)2 8509 5401  
Fax: +61 (0)2 8509 5403  
Email: [Australiansales@bestglove.com](mailto:Australiansales@bestglove.com)