

# SAFEGRIP®

POWDER-FREE LATEX EXAMINATION GLOVES

• FOR HIGH-RISK ENVIRONMENTS



## *Recommended for use when:*

- Added protection is required for extreme environments
- Extra coverage over the wrist and forearm is needed for enhanced safety and security
- Working in situations requiring intense pressure and strain on the hands



For the highest level of protection and durability in latex exam gloves for heavy-duty or high-risk environments, rely on **SafeGrip®** gloves.

**SafeGrip®** gloves are up to 4 times thicker than standard latex exam gloves for defense against many hazards, and is ideal protection for your toughest tasks.

**SafeGrip®** gloves are the best choice in a latex for heavy-duty/high-risk applications.

Glove meets or exceeds the emergency medical glove requirements of NFPA 1999, *Standard on Protective Clothing for Emergency Medical Operations*, 2003 edition.\*

\*Certification on file.

**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®

REGISTERED  
5013485

# SAFEGRIP®

## POWDER-FREE LATEX EXAMINATION GLOVES

### ● FOR HIGH-RISK ENVIRONMENTS



#### GLOVE PROFILE

	ASTM*	SafeGrip®
<b>Length</b> (mm/in):	220/8.7	300/11.8
<b>Cuff Thickness</b> (mils):	n/a	7.9
<b>Palm Thickness</b> (mils):	min. 3.1	11.4
<b>Finger Thickness</b> (mils):	min. 3.1	16.5
<b>Tensile Strength</b>		
Before Aging (MPa):	min. 18	21
After Aging (MPa):	min. 14	18
<b>Elasticity/Elongation:</b>		
Before Aging (%):	min. 650	750
After Aging (%):	min. 500	600
<b>Reduced Protein Rating:</b>	max. 200 µg/dm <sup>2</sup>	max. 50 µg/g

#### GLOVE ATTRIBUTES

<b>Material:</b>	Premium quality natural rubber latex
<b>Type:</b>	Non-sterile
<b>Shape:</b>	Ambidextrous
<b>Color:</b>	Blue
<b>Sizes:</b>	Small, Medium, Large, X-Large
<b>Cuff:</b>	Extended, beaded
<b>Glove Interior:</b>	Powder-free
<b>Glove Exterior:</b>	Textured fingertips
<b>Usage:</b>	For single use only
<b>Packaging:</b>	50 gloves per dispenser box, by weight 10 boxes per case

\* ASTM D 3578

See for yourself why **Microflex®** gloves are the first choice of exam glove users nationwide. For samples or more information on any of the fine latex and synthetic products from **Microflex®**, contact your **Microflex®** distributor or for the name of a distributor in your area, call:

**(800) 876-6866**

Or visit our website at:

**[www.microflex.com](http://www.microflex.com)**

Fax: (800) 876-6632.

PRODUCT CODE	RANGE OF SIZE	DESCRIPTION	CASE QUANTITY
SG-375	SMALL	Latex, Non-sterile, Examination	500
SG-375	MEDIUM	Latex, Non-sterile, Examination	500
SG-375	LARGE	Latex, Non-sterile, Examination	500
SG-375	EXTRA-LARGE	Latex, Non-sterile, Examination	500



**CAUTION STATEMENT:** This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.



**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®

2301 ROBB DRIVE • RENO, NV 89523