

SemperCARE®

by Sempermed®

NITRILE EXAMINATION GLOVES

Industrial

Exam

Surgical

Non-Latex • Powder-Free • Textured



Sempermed®

We are people protection.™

Sempermed is a vertically integrated manufacturer that provides focused sales, marketing, and logistical support to those in need of hand protection.



SemperCare® Nitrile is manufactured as a copolymer glove to provide a better fit and feel. This non-latex alternative is made stronger through a manufacturing process called “fusion bonding” – similar to our Polymed® glove. “Fusion bonding” integrates synthetic polymers into the material to create a multi-layered glove for enhanced protection.

FEATURES

- “Fusion bonding”
- Better elasticity
- Powder-Free
- Tested for protection against chemotherapy drugs
- Manufactured without natural rubber

BENEFITS

- Strong and more durable
- Less pressure on fingers for an improved fit with more comfort
- Eliminates aerosolization of glove powder
- Appropriate for handling cytotoxic chemicals
- Eliminates exposure to latex

GLOVE COUNT (By Weight)

200 Gloves/Box (180 Gloves/Box XL)

10 Boxes/Case

THICKNESS* (Single Wall)

Finger: 0.05 mm/2.0 mils min.

Palm: 0.05 mm/2.0 mils min.

PHYSICAL PROPERTIES*

Before aging:

Tensile Strength: 14 MPa min.

Elongation: 500% min.

*The specifications provided reflect the minimum acceptable requirements of ASTM and FDA where applicable. Gloves offered by Sempermed USA, Inc. meet or exceed these physical requirements. Independent laboratory test results are available upon request. All Sempermed gloves are single use only.

PRODUCT CHARACTERISTICS

- Synthetic Rubber Copolymer Nitrile Blend
- Textured Fingertips
- Powder-Free
- Ambidextrous
- Beaded Cuff
- Single Use Only

Gloves used for protection against chemotherapy drug exposure must be selected specifically for the type of chemicals used. Review material safety data sheets for the drug being used to determine the required level of protection. Test results available upon request. Chemotherapy drugs tested.

| Order # | Size | Length (min.) | Palm Width (min.) |
|----------|------|---------------|-------------------|
| NIPFT201 | XS | 220mm | 70mm +/- 10mm |
| NIPFT202 | S | 220mm | 80mm +/- 10mm |
| NIPFT203 | M | 230mm | 95mm +/- 10mm |
| NIPFT204 | L | 230mm | 110mm +/- 10mm |
| NIPFT205 | XL | 230mm | 120mm +/- 10mm |

Sempermed®

We are people protection.™

www.SempermedUSA.com • 13900 49th Street North • Clearwater, FL 33762

Phone: 1.800.366.9545 • Fax: 1.800.763.5491

Outside the U.S.: 727.787.7250 • Fax: 727.787.7558

