

SemperCARE®

by Sempermed®

PVC EXAMINATION GLOVES

Industrial

Exam

Surgical

Non-Latex • Powder-Free



Sempermed®

We are people protection.™

Sempermed is a vertically integrated manufacturer that provides focused sales, marketing, and logistical support to those in need of hand protection.



Strength. Durability. Sensitivity. Fit. The SemperCare™ PVC glove is designed as the ultimate solution for those who are at risk of becoming sensitive to natural rubber latex. It's cream color and new resin is an upgrade from regular vinyl gloves, giving the SemperCare™ glove the look and feel of latex without the latex allergens. It offers exceptional strength and protection, yet maintains maximum tactile sensitivity.

FEATURES

- Manufactured without natural rubber
- Powder-Free
- Cream-color vinyl
- Supple touch

BENEFITS

- A latex alternative
- Eliminates glove powder aerosolization
- Latex glove appearance and feel
- Maximum tactile sensitivity

GLOVE COUNT (By Weight)

100 Gloves/Box (90 Gloves/Box XL)
10 Boxes/Case

THICKNESS* (Single Wall)

Finger: 0.05 mm/2.0mils min.
Palm: 0.08 mm/3.2mils min.

PHYSICAL PROPERTIES*

Before aging:
Tensile Strength: 11 MPa min.
Elongation: 300% min.

*The specifications provided reflect the minimum acceptable requirements of ASTM and FDA where applicable. Gloves offered by Sempermed USA, Inc. meet or exceed these physical requirements. Independent laboratory test results are available upon request.
All Sempermed gloves are single use only.

PRODUCT CHARACTERISTICS

- Polyvinyl Chloride Blend (PVC)
- Powder-Free
- Ambidextrous
- Beaded Cuff
- Single Use Only

Order #	Size	Length (min.)	Palm Width (avg.)
SCVNP101	XS	230mm	75mm +/- 5mm
SCVNP102	S	235mm	85mm +/- 5mm
SCVNP103	M	235mm	95mm +/- 5mm
SCVNP104	L	235mm	105mm +/- 5mm
SCVNP105	XL	235mm	115mm +/- 5mm

Sempermed®

We are people protection.™

www.SempermedUSA.com • 13900 49th Street North • Clearwater, FL 33762
Phone: 1.800.366.9545 • Fax: 1.800.763.5491
Outside the U.S.: 727.787.7250 • Fax: 727.787.7558

