

SEMPERFORCE™

BLACK NITRILE EXAMINATION GLOVES

Industrial

Exam

Surgical


Non-Latex • Powder-Free • Textured



Sempermed® 

We are people protection.™

Sempermed is a vertically integrated manufacturer that provides focused sales, marketing, and logistical support to those in need of hand protection.



SemperForce™ black nitrile examination gloves by Sempermed are designed for demanding work environments where your safety and comfort are a must. A textured surface promises secure gripping capabilities. SemperForce™ black nitrile gloves are ideal for automotive, industrial, tattoo, and cosmetology applications.

FEATURES

- Black color
- Manufactured without natural rubber
- Textured fingertips

BENEFITS

- Hides stains caused from grease, dirt, dyes, and ink
- Eliminates exposure to latex and potential allergens
- Enhanced dexterity and gripping capabilities

GLOVE COUNT (By Weight)

100 Gloves/Box (90 Gloves/Box XL-XXL)
10 Boxes/Case

THICKNESS* (Single Wall)

Finger: 0.05 mm/2.0 mils min.
Palm: 0.05 mm/2.0 mils min.

PHYSICAL PROPERTIES*

Before Aging:

Tensile Strength: 14 MPa min.
Elongation: 500% min.

*The specifications provided reflect the minimum acceptable requirements of ASTM and FDA where applicable. Gloves offered by Sempermed USA, Inc. meet or exceed these physical requirements. Independent laboratory test results are available upon request.
All Sempermed gloves are single use only.

PRODUCT CHARACTERISTICS

- Synthetic Nitrile Rubber (NBR)
- Textured Fingertips
- Powder-Free
- Ambidextrous
- Beaded Cuff
- Single Use Only

| Order # | Size | Length (min.) | Palm Width (min.) |
|---------|------|---------------|-------------------|
| BKNF102 | S | 230mm | 85mm +/- 5mm |
| BKNF103 | M | 230mm | 95mm +/- 5mm |
| BKNF104 | L | 230mm | 108mm +/- 5mm |
| BKNF105 | XL | 230mm | 113mm +/- 5mm |
| BKNF106 | XXL | 230mm | 118mm +/- 5mm |

Sempermed®
We are people protection.™

www.SempermedUSA.com • 13900 49th Street North • Clearwater, FL 33762
Phone: 1.800.366.9545 • Fax: 1.800.763.5491
Outside the U.S.: 727.787.7250 • Fax: 727.787.7558