

Synthetic

by **Sempermed**®

VINYL EXAMINATION GLOVES

Industrial

Exam

Surgical

Non-Latex • Powder-Free • Smooth



Sempermed®

We are people protection.™

Sempermed is a vertically integrated manufacturer that provides focused sales, marketing, and logistical support to those in need of hand protection.



Synthetic brand examination gloves by Sempermed® provide a solution for professionals and their patients who are at risk of becoming sensitive to natural rubber latex. They maximize tactile sensitivity and provide a contoured fit, reducing hand fatigue.

FEATURES

- Manufactured without natural rubber
- Powder-free
- Supple touch

BENEFITS

- Eliminates exposure to latex
- Eliminates aerosolization of glove powder
- Maximum tactile sensitivity

GLOVE COUNT (By Weight)

100 Gloves/Box (90 Gloves/Box XL)
10 Boxes/Case

THICKNESS* (Single Wall)

Finger: 0.05 mm/2.0 mils min.
Palm: 0.08 mm/3.2 mils min.

PHYSICAL PROPERTIES*

Before aging:

Tensile Strength: 11 MPa min.

Elongation: 300% min.

*The specifications provided reflect the minimum acceptable requirements of ASTM and FDA where applicable. Gloves offered by Sempermed USA, Inc. meet or exceed these physical requirements. Independent laboratory test results are available upon request.

All Sempermed gloves are single use only.

PRODUCT CHARACTERISTICS

- Polyvinyl Chloride (PVC)
- Powder-Free
- Ambidextrous
- Single Use Only

Order #	Size	Length (min.)	Palm Width (avg.)
EVNP101	XS	230mm	76mm +/- 5mm
EVNP102	S	230mm	85mm +/- 5mm
EVNP103	M	230mm	95mm +/- 5mm
EVNP104	L	230mm	105mm +/- 5mm
EVNP105	XL	230mm	115mm +/- 5mm

Sempermed®

We are people protection.™

www.SempermedUSA.com • 13900 49th Street North • Clearwater, FL 33762

Phone: 1.800.366.9545 • Fax: 1.800.763.5491

Outside the U.S.: 727.787.7250 • Fax: 727.787.7558