


SUPRENO[®] EC

POWDER-FREE NITRILE EXAM GLOVES

• FOR HIGH RISK ENVIRONMENTS

Recommended for use when:

- Strength and protection are of great importance
- Day-to-day hazards require the best combination of durability and security
- Working in environments that demand consistent grip and reliability



Exceptional Protection in a Tough Nitrile Glove.

Supreno[®] EC is a durable nitrile glove for high-risk applications. **Supreno[®] EC** gloves offer all the benefits of a nitrile glove, with the most important being unique and reliable protection against the hazards you work with. **Supreno[®] EC** gloves are formulated for enhanced strength and durability while maintaining a comfortable fit.

Supreno[®] EC gloves meet or exceed the emergency medical glove requirements of NFPA 1999, *Standard on Protective Clothing for Emergency Medical Operations*, 2003 edition.*

*Certification on file.

NOT MADE FROM NATURAL RUBBER LATEX

MICROFLEX[®]
THE MOST TRUSTED NAME IN GLOVES[®]

REGISTERED
ISO 13485

SUPRENO[®] EC

POWDER-FREE NITRILE EXAM GLOVES

• FOR HIGH RISK ENVIRONMENTS

GLOVE PROFILE

	ASTM*	Supreno [®] EC
Length (mm/in):	220/8.7	295/11.6
Cuff Thickness (mils):	n/a	4.7
Palm Thickness (mils):	min. 2.0	5.9
Finger Thickness (mils):	min. 2.0	7.9
Tensile Strength		
Before Aging (MPa):	min. 14	min. 14
After Aging (MPa):	min. 14	min. 14
Elasticity/Elongation:		
Before Aging (%):	min. 500	min. 500
After Aging (%):	min. 400	min. 400

GLOVE ATTRIBUTES

Material:	Nitrile; not made from natural rubber latex
Type:	Non-sterile
Shape:	Ambidextrous
Color:	Blue
Sizes:	Small, Medium, Large, XL, XXL, XXXL
Cuff:	Extended, beaded
Glove Interior:	Polymer-coated, powder-free
Glove Exterior:	Textured fingers
Usage:	For single use only
Packaging:	S-XXL: 50 gloves per dispenser box, by weight XXXL: 40 gloves per dispenser box, by weight 10 boxes per case

* ASTM D 6319



Nitrile Examination Gloves

Exceptional Comfort in High Risk Hand Protection.

See for yourself why **Microflex[®]** gloves are the first choice of exam glove users nationwide. Contact a **Microflex[®]** distributor today to purchase quality **Microflex[®]** products. For the name of a distributor in your area, please call **(800) 876.6866**.

Or visit our website at:

www.microflex.com

Fax: (800) 876-6632.



PRODUCT CODE	RANGE OF SIZE	DESCRIPTION	CASE QUANTITY
SEC-375-S	SMALL	Nitrile, Non-sterile, Examination	500
SEC-375-M	MEDIUM	Nitrile, Non-sterile, Examination	500
SEC-375-L	LARGE	Nitrile, Non-sterile, Examination	500
SEC-375-XL	EXTRA-LARGE	Nitrile, Non-sterile, Examination	500
SEC-375-XXL	DOUBLE EXTRA-LARGE	Nitrile, Non-sterile, Examination	500
SEC-375-3XL	TRIPLE EXTRA-LARGE	Nitrile, Non-sterile, Examination	400

Caution Statement: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.

U.S. PATENT NO. RE. 35.616. EUROPEAN PATENT NO. 0456333.

©2007 MICROFLEX CORP. ALL RIGHTS RESERVED. 07.05536C MRK AL 101707 AL

MICROFLEX[®]
THE MOST TRUSTED NAME IN GLOVES[®]
2301 ROBB DRIVE • RENO, NV 89523