

Microflex®  
**TRANQUILITY®**

**Powder-Free White  
 Nitrile\* Exam Gloves**

**Item #**

- Extra Small..... TQ-601-XS
- Small ..... TQ-601-S
- Medium..... TQ-601-M
- Large..... TQ-601-L
- Extra Large..... TQ-601-XL



**Recommend for use in the following industries:**  
 Healthcare / Dental / Industrial Safety / Laboratory

**Rejuvenate your hands with Tranquility® gloves!**

For the perfect balance of comfort and protection, choose **Tranquility®** white nitrile. **Tranquility®** gloves are latex-free and contain no natural rubber proteins. The glove of choice for individuals that require a soft, thin, comfortable glove for everyday use. **Tranquility®** gloves conform to the hand allowing for excellent tactile sensitivity and grip. **Tranquility's** compact box size saves on storage requirements.

**Recommended for use when:**

- A latex-free solution is required, at a value price
- Enhanced comfort and tactile sensitivity are a must
- Everyday use of small instruments demands superior wet and dry grip

**NOT MADE FROM  
 NATURAL RUBBER LATEX**

\*U.S. Patent No. 6,451,893, U.S. Patent No. 7,176,260

Microflex® Tranquility® Profile			For current specification information, please visit <a href="http://www.microflex.com">www.microflex.com</a> .	
	ASTM D 6319	Tranquility®		
<b>Length (mm/in):</b>	220 / 8.7	240 / 9.6	<b>Material:</b> Nitrile	<b>Type:</b> Non-sterile
<b>Cuff Thickness (mils):</b>	n/a	2.8 ± 0.8	<b>Size:</b> XS - XL	<b>Color:</b> White
<b>Palm Thickness (mils):</b>	min. 2.0	3.1 ± 0.5	<b>Shape:</b> Ambidextrous	<b>Usage:</b> For single use only
<b>Finger Thickness (mils):</b>	min. 2.0	4.3 ± 1.2	<b>Cuff:</b> Standard examination length (9.5"), Beaded	
<b>Tensile Strength:</b>			<b>Glove Exterior:</b> Textured fingers	
Before Aging (MPa):	min. 14	min. 18	<b>Glove Interior:</b> Powder-free, chlorinated	
After Aging (MPa):	min. 14	min. 16	<b>Packaging:</b> 1,000 gloves/case, by weight; 10 boxes/case	
<b>Elasticity/Elongation:</b>			100 gloves per dispenser box, by weight	
Before Aging (%):	min. 500	min. 500	<b>Caution:</b> Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.	
After Aging (%):	min. 400	min. 400		

**The reason is the quality... the difference is Microflex®.**

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit [www.microflex.com](http://www.microflex.com). For the name of distributor in your area, please **CALL (1-800-876-6866)** 6am-5pm (pst) or **FAX (1-800-876-6632)** 24 hours

