


Microflex® XCEED™

Powder-Free Nitrile Examination Gloves

Item

Extra-Small.....	XC-310-XS
Small	XC-310-S
Medium	XC-310-M
Large	XC-310-L
Extra-Large	XC-310-XL*



Recommended for use in the following industries:

Dental / Healthcare / Laboratory / Tattoo / Industrial Safety / Manufacturing / Automotive

New Revolutionary Nitrile Formulation and Manufacturing Process delivers enhanced comfort, strength and value.

Exceeding your expectations...

Comfort... XCEED™ gloves are soft and provide improved tactile sensitivity which helps reduce hand fatigue.

Strength... New formulation raises the bar in strength and durability in disposable examination gloves.

Value Packed... Each box contains 250 gloves which is 2.5 times more gloves than the average 100 count box. Value packing saves you money and time. In addition, value packing is better for the environment.

- The comfort of latex without the risk of latex sensitivities
- Enhanced tactile sensitivity and reduced hand fatigue
- Raises the bar in strength and durability
- Value packed 250* gloves per box

Try XCEED™ gloves today and exceed your own expectations in comfort, strength and value.

www.microflex.com/XCEED

Microflex® XCEED™ Profile

For current specification information, please visit www.microflex.com.

	ASTM D 6319	XCEED™		
Length (mm/in):	220 / 8.7	235 / 9.3	Material: Nitrile (Latex-Free)	Type: Non-sterile
Cuff Thickness (mils):	n/a	2.4	Size: XS - XL	Color: Blue
Palm Thickness (mils):	min. 2.0	2.8	Shape: Ambidextrous	Usage: For single use only
Finger Thickness (mils):	min. 2.0	3.5	Cuff: Beaded	
Tensile Strength:			Glove Exterior: Textured fingertips	
Before Aging (MPa):	min. 14	22	Glove Interior: Powder-free	
After Aging (MPa):	min. 14	14	Packaging: 2,500 gloves/case, by weight; 10 boxes/case	
Elasticity/Elongation:			250 gloves per dispenser box, by weight	
Before Aging (%):	min. 500	500	*Extra Large: 230 gloves per dispenser box, by weight	
After Aging (%):	min. 400	500	<small>Caution: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.</small>	

The reason is the quality.. the difference is Microflex®.

See for yourself why Microflex® gloves are the first choice for exam glove users nationwide. For more information, or to order a sample, please call: (800) 876-6866 or visit: www.microflex.com/XCEED

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
A subsidiary of BarrierSafe Solutions International, Inc.

+ [STRENGTH]

NEW! Revolutionary Nitrile Formulation delivers up to 60% stronger material than leading brands.

+ [COMFORT]

NEW! Manufacturing Processes offers you a soft and comfortable nitrile glove that rivals the fit of latex.

ENGINEERED PERFORMANCE

Are XCEED™ gloves the best gloves yet? *We think so!*

Our skillfully engineered **NEW** polymer technology delivers the right comfort, the right strength all at the right price.

+ [VALUE]

NEW! Value Packing enables XCEED™ gloves to contain 250* gloves per box and save you money, time, storage space and is better for the environment.

NEW!



MICROFLEX
THE MOST TRUSTED NAME IN GLOVES®

a subsidiary of BarrierSafe Solutions International, Inc.