



Best® Latex 5005(PF)

**100% NATURAL RUBBER LATEX
DISPOSABLE GLOVE
5 MILS, 9.5 INCHES
LABORATORY GLOVE**

Description

SHOWA Best® Latex 5005PF Disposable Natural Rubber Latex glove is specially suited to customers looking for a low cost, high quality disposable Latex glove.

Best® Latex 5005PF (powder-free) is a five-mil, 9.5-inch disposable glove that is a perfect solution for customers who require the barrier protection found in latex gloves. The glove is ideally suited for many laboratory or industrial applications such as automotive, janitorial/maintenance and food service.



Textured Fingertips for enhanced grip

Powder-Free Disposable Latex features Low Modulus formulation for Maximum Dexterity

Protection from aqueous and powdered Laboratory and Industrial Chemicals

Available in a wide range of sizes for every hand

Sizes:

- 5005PFS small**
- 5005PFM medium**
- 5005PFL large**
- 5005PFXL X-large**

Suggested Applications:

- Janitorial
- General Maintenance
- Assembly
- Housekeeping
- Laboratories
- Food Service
- Parts Assembly

This list of suggested applications is partial and based on laboratory testing of glove performance. For a complete list of suggested applications, visit our website: www.showabestglove.com



Packaging:

100 gloves per dispenser
10 dispensers per case

in compliance 21 CFR, parts 170 to 199